50P13 Manu Senza

DE	Gebrauchsanweisung	3
EN	Instructions for Use	7
FR	Instructions d'utilisation	12
IT	Istruzioni per l'uso	17
ES	Instrucciones de uso	22
PT	Manual de utilização	27
NL	Gebruiksaanwijzing	31
SV	Bruksanvisning	36
DA	Brugsanvisning	41
NO	Bruksanvisning	45
FI	Käyttöohje	50
PL	Instrukcja użytkowania	54
HU	Használati utasítás	59
CS	Návod k použití	64
SK	Návod na používanie	68
TR	Kullanma talimatı	73
EL	Οδηγίες χρήσης	77
RU	Руководство по применению	82
CN	使用说明书	87

1

2

3

4

	Side / Seite	Size / Größe	Wrist circumference / Handgelenkumfang	
	L/R	XS	14 – 15 cm	5.5 – 5.9 inch
	L/R	S	15 – 16 cm	5.9 – 6.3 inch
	L/R	M	16 – 17 cm	6.3 – 6.7 inch
	L/R	L	17 – 18 cm	6.7 – 7.1 inch
	L/R	XL	18 – 19 cm	7.1 – 7.5 inch
	L/R	XXL	19 – 20 cm	7.5 – 7.9 inch

Art.-Nr.: 50P13=R(Side)-M(Size)-1

Material	Knit Fabric / Gestrick: PA, EL
	Terry Cloth Covers / Frotteebezüge: Cotton / Baumwolle, PES
	Plastic Splint / Kunststoffschiene: PA
	Wrist Strap / Handgelenkgurt: PA, EL
	Pad / Pelotte: viscoelastic Technogel® / viskoelastisches Technogel®

1 Verwendung

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1 Application

1.1 Intended Use

The support / orthosis is to be used **exclusively** for orthotic wrist fittings and is **only** suitable for contact with intact skin.

1.2 Effects

The form-knitted support / orthosis provides targeted wrist compression and stabilisation. The support / orthosis promotes sensorimotor functions and activates muscular stabilisation. The anatomically pre-shaped and individually adaptable volar plastic splint keeps the hand in an ideal position and limits wrist movement. The pad in the radial section of the wrist enhances the micro-massage and promotes the resorption of oedemas, haematomas and effusions.

1.3 Indications

- Irritation (e. g. tendovaginitis, tendomyopathy, minor joint effusion, post-traumatic [distorsion], post-operative, reactive irritation and arthralgia)
- Rheumatoid arthritis
- Carpal tunnel syndrome
- Degenerative wrist disorders (arthrosis)
- Feelings of instability (e. g. snapping ulna)

The indications must be determined by the physician.

1.4 Contraindications

1.4.1 Absolute Contraindications

None.

1.4.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases / injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb; lymphatic flow disorders, including

unclear soft tissue swelling distant to the body area to which the support / orthosis will be applied; sensory and circulatory disorders in the body region being treated and in the area of the fingers.

1.5 Additional Restrictions for Use

The support / orthosis is designed for use on **one** patient. Re-use is not permitted. The support / orthosis must be used in keeping with the indications and according to the instructions of the physician.

2 Safety

2.1 Explanation of warning symbols

⚠ WARNING	Warnings regarding possible risks of severe accident or injury.
⚠ CAUTION	Warnings regarding possible risks of accident or injury.
NOTICE	Warnings regarding possible technical damage.

2.2 General Safety Instructions

⚠ CAUTION

The support / orthosis may be over-tightened

Risk of local pressure and constriction of blood vessels and nerves.

- ▶ Verify the correct positioning and fit of the support / orthosis.

⚠ CAUTION

Use on more than one patient

Skin irritation, formation of eczema or infections

- ▶ Only use the support / orthosis on one patient.

NOTICE

Incorrect environmental conditions

The support / orthosis is not flameproof.

- Do not expose the support / orthosis to open flame or other heat sources.

NOTICE

The lifespan of the materials may be reduced

- Avoid contact with substances containing oils and acids, salves and lotions.

INFORMATION

Improper changes to the support / orthosis are not permitted and will eliminate product liability.

2.3 Information to be Provided to the Patient

NOTICE

Damage due to wear and tear

- Instruct the patient to inspect the support / orthosis for proper functionality, wear and damage prior to each use.
- Advise the patient that the support / orthosis must not be used if a closure or other part of the support / orthosis shows signs of wear.

INFORMATION

The patient must be instructed to see a physician immediately if any exceptional changes are noted (e. g. worsening of the complaint).

3 Handling

INFORMATION

The patient must be instructed in the care and handling of the support / orthosis. The initial fitting and application of the support / orthosis must be carried out by trained, qualified personnel. The daily duration of use and the total period of application are generally determined by the physician.

3.1 Size Selection

- 1) Determine the wrist circumference.
- 2) Select the size of the support (see sizing table).

3.2 Adaptation and Application

The volar plastic splint is anatomically pre-shaped. In the course of the fitting, the plastic splint can be heated and then individually shaped by the orthotist. The support is delivered closed and ready for trial fitting.

- 1) Open the elastic strap.
- 2) Insert the plastic splint into the pocket on the inside of the support (see Fig. 1). Ensure that the lower end of the plastic splint is in the pocket intended for this purpose.
- 3) Pull the support over the hand (see Fig. 2), taking into account the opening provided for the thumb (see Fig. 3).
- 4) Close the elastic strap (see Fig. 4).
- 5) The elastic strap can be removed if required.
- 6) When providing the support to the patient, the proper fit and positioning of the support must be verified.

3.3 Cleaning

Clean the support / orthosis regularly.

- 1) Remove the plastic splint.
- 2) **Recommendation:** put the support / orthosis in a bag or net for washing.

- 3) Wash the support / orthosis in 30 °C warm water using a standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e. g. sunlight, stove or radiator).

4 Legal information

4.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

4.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH