
50R59 Dorso Direxa Posture


DE Gebrauchsanweisung	4
EN Instructions for Use	9
FR Instructions d'utilisation	14
IT Istruzioni per l'uso	19
ES Instrucciones de uso	24
NL Gebruiksaanwijzing	30
PL Instrukcja użytkowania	35
HU Használati utasítás	41
CS Návod k použití	46
SK Návod na používanie	51
TR Kullanma talimatı	56
RU Руководство по применению	61
CN 使用说明书	66

 Size / Größe	Waist circumference / Tailenumfang*	
	cm	inch
XS	70 – 80	27.6 – 31.5
S	80 – 90	31.5 – 35.4
M	90 – 100	35.4 – 39.4
L	100 – 110	39.4 – 43.3
XL	110 – 120	43.3 – 47.2
XXL	120 – 130	47.2 – 51.2

*approx. 2 cm/0.8 inch above the iliac crest / ca. 2 cm oberhalb des Beckenkamms

Material	Fabric / Gewebe: PA, EL, Lycra®
	Terry Cloth Covers / Frotteebezüge: Cotton / Baumwolle, PES
	Corset Splints / Korsettstäbe: Metal / Metall, POM
	Pad / Pelotte: TPE

1 Verwendung

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1 Application

1.1 Intended Use

The orthosis is intended **exclusively** for orthotic fittings of the trunk and **exclusively** for contact with intact skin.

1.2 Effects

The orthosis activates muscular stabilisation of the thoracic spine and straightens it actively. It promotes sensorimotor functions, supports pain relief, and relieves the thoracic and lumbar spine.

1.3 Indications

- Poor body posture (hunchback without fused vertebrae, Scheuer-mann's disease)
- Osteoporosis of the thoracic and lumbar spine
- Aftercare of osteoporotic fractures of the thoracic and/or lumbar spine
- Back pain (thoracic and lumbar spine)
- Slight degenerative changes of the thoracic and lumbar spine (e.g. osteochondrosis, facet joint syndrome)

Indications must be determined by the physician.

1.4 Contraindications

1.4.1 Absolute Contraindications

None.

1.4.2 Relative contraindications

The following indications require consultation with a physician: Skin diseases/injuries, inflammatory changes, raised scars with swelling, reddening and warming of the treated body areas, tumours, lymphatic flow disorders - including unclear soft tissue swelling in areas distant to where the device is applied, sensitivity disorders in the arms and legs.

1.5 Further usage restrictions

The orthosis is designed for use on **one** patient. Re-use is not permissible. The orthosis must be used according to the indications and following the instructions of a physician.

2 Safety

2.1 Explanation of warning symbols

⚠ WARNING	Warnings regarding possible risks of severe accident or injury.
⚠ CAUTION	Warnings regarding possible risks of accident or injury.
NOTICE	Warnings regarding possible technical damage.

2.2 General Safety Instructions

⚠ CAUTION

The support / orthosis may be over-tightened

Risk of local pressure and constriction of blood vessels and nerves.

- ▶ Verify the correct positioning and fit of the support / orthosis.

⚠ CAUTION

Use on more than one patient

Skin irritation, formation of eczema or infections

- ▶ Only use the support / orthosis on one patient.

NOTICE

Incorrect environmental conditions

The support / orthosis is not flameproof.

- Do not expose the support / orthosis to open flame or other heat sources.

NOTICE

The lifespan of the materials may be reduced

- Avoid contact with substances containing oils and acids, salves and lotions.

INFORMATION

Improper changes to the support / orthosis are not permitted and will eliminate product liability.

2.3 Notes on handover to the patient

△ CAUTION

Direct skin contact with the support/orthosis can cause skin irritation

- Inform the patient that he/she should not wear the support/orthosis directly against the skin.

NOTICE

Damage due to wear and tear

- Instruct the patient to inspect the support / orthosis for proper functionality, wear and damage prior to each use.
- Advise the patient that the support / orthosis must not be used if a closure or other part of the support / orthosis shows signs of wear.

INFORMATION

The patient must be instructed to see a physician immediately if any exceptional changes are noted (e. g. worsening of the complaint).

3 Handling

INFORMATION

The daily duration of use and usage period are generally determined by the physician. Initial fitting and application of the support/orthosis must be performed by trained specialists. The patient must be instructed on how to handle and care for the support/orthosis.

3.1 Selecting the size

- 1) Measure the waist circumference approx. 2 cm above the iliac crest.
- 2) Determine the orthosis size (see size table).

3.2 Fitting and Application

- 1) **Fitting the metal rods** (see Fig. 1): Remove the 2 metal rods from the pockets and fit them to the spine individually. Put the metal rods back in the pockets.
- 2) **Application:** Place the upper part of orthosis over the shoulders like a backpack (see Fig. 2). Place the lower part of the orthosis around the trunk. Place the fingers in the grip pockets of the abdominal closure. Apply equal force to both sides to pull the support forward (see Fig. 3).
- 3) **Hook and loop material:** Attach the right part to the top of the left part of the abdominal closure (see Fig. 4). Remove the fingers from the left and then from the right grip pockets. Ensure that the hook and loop material sections lie flat against each other and are not in contact with clothing or skin. Route the ends of the shoulder strap to the front around the chest (see Fig. 5) and secure the hook and loop closure below the chest (see Fig. 6)
- 4) **Increasing compression:** Pull the ends of the elastic straps forward (see Fig. 7) and attach them to the bottom closure (see Fig. 8).
- 5) **Optional – install the pad** (see Fig. 9): Use the hook and loop closure to attach the pad to the inside of the orthosis so that it is at least 1 cm above the seat surface when seated.

- 6) **Final check:** The correct fit and proper positioning of the orthosis must be checked during handover to the patient. The orthosis is seated correctly when the lower edge of the orthosis extends just below the groin crease.
- 7) **Instruct the patient:** Explain to the patient how to don the orthosis and ensure that it is seated correctly. Practice donning and doffing of the orthosis with the patient.

3.3 Cleaning

Clean the orthosis regularly.

- 1) Remove the metal rods. Close all hook and loop closures.
- 2) **Recommendation:** Use a laundry bag or net.
- 3) Wash the orthosis in 40 °C warm water with standard mild detergent.
Do not use fabric softener.
- 4) Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).

4 Disposal

Dispose of the product in accordance with national regulations.

5 Legal information

5.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

5.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH