Generation Osteoporose-Orthese

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Relief and support is a basic principle when treating injuries and excess stress as they considerably accelerate the healing process and prevent resulting health damages.


The BORT Generation Osteoporosis Brace counteracts kyphosis in the thoracic spine region by activating back muscle, improves the posture and supports the strengthening of the muscles.

Instructions for use for the patient:

The brace should be worn directly over the clothes. Pay attention to avoid creasing the clothes in the vicinity of the brace as this can result in pressure sores.

1. Opening the front, wide brace fastener.
2. Guide the arms through the shoulder loops as accustomed for a rucksack.
3. Then close the front, wide fastener on the stomach under slight tension. Ensure a good fit, the ventral plate should be positioned a little above the groin.
4. Straighten yourself up actively ("make yourself big") and tighten both shoulder belts in sync by pulling them lightly.
5. If you find the correction unpleasant, reduce it by tilting the shoulder belts' guide loops downwards and upwards with your thumbs as for a rucksack.
6. To remove the brace, release the shoulder belts as described above, open the front stomach fastener and take the brace off like a rucksack.

- Information:**
- Should you incur pressure sores, pain or other problems whilst wearing the brace, remove it and contact your technician.
 - Allow yourself time to become accustomed to the brace and increase the wearing time gradually.
 - If you have questions regarding the use and wearing time, please refer to your doctor or technician.

Indications:

Active straightening for relief and correction of the lumbar / thoracic spine associated with restricted movement in the sagittal plane, i.e. stable, osteoporotic vertebral fractures, osteoporosis, cautionary treatment of secondary kyphosis (stable fractures, tumours), muscular insufficiency, pre-operative and post-operative stabilization

Contraindications:

Inflammatory or injury-caused skin changes (including swellings and redness) in the body region to be treated. Circulatory problems, soft tissue swellings. Neurologically caused sensory disorders and skin trophism in the body region treated (sensitivity disorders, sensory disorders with and without skin damages). Ongoing immobilisation (contracture), in particular with older people

Composition:

Base: 42 % polyurethane, 35 % polyamide, 10 % cotton,
7 % polyester, 5 % elastane, 1 % Helanca

Belts/fasteners: 60 % polyamide, 40 % polyurethane

Important Note:

No re-use – this medical aid is intended for treating one patient. If wearing at night is prescribed, avoid negative influence on the circulatory system. In case of numbness, loosen the brace immediately and remove it if necessary. In case of persistent complaints, consult the doctor. Do not use any cream or ointment when wearing the brace as they can damage the material.


Information for the technician:

Information: The brace may only be adapted by expert staff!


Adaption of the brace on the patient:

1. Shape the aluminium back splint to the contours of the patient's spine. Ensure sufficient distance of the back plate to the body on the upper (cranial) edge to the body (approx. 2 cm). Only with this distance is it possible to straighten the patient. If required, the back splint can be removed from the textile cover.

2. Place the brace on the patient, it can be positioned by the patient like a rucksack, then close the ventral plate. Ensure correct positioning of the brace, the lower edge must be at the height of the gluteal fold. The front ventral plate should be positioned at the height of the pubic bone and not restrict hip flexion.


3. Then set the brace belt lengths. Remove the crocodile clasp ends of the belts and shorten the strap lengths accordingly. Start with the belt straps on the hip, always tighten the belts on both sides in sync. After completion, tighten the crocodile clasp ends and close the belt strap under slight tension. Pay attention to an appropriate length.

4. Next, adjust the combined torso and shoulder belts and proceed as described above. Before finally shortening the belts, ensure that there is a gap of approx. 3 cm between the pull loops and the guide loops for the shoulder belts.


5. Pad the shoulder belts in the armpit region using the enclosed belt pads.

6. To reduce tension to the armpits and avoid strangulation, you can correct the course of the belt by changing the position of the dorsal guide loops (attached to the back plate with Velcro). The level of the costal arch serves as a reference point.

7. In order to achieve ideal straightening of the patient, then use the pull loops in sync to tighten the torso and shoulder belts. Avoid overcorrection!

To reduce tension or remove the brace, tilt the shoulder belts' guide loops downwards and upwards with your thumbs as for a rucksack.

8. Check the correct fit of the brace.

- Are the belt straps too tight or do they cut in?
- Is adequate space (approx. 2 cm) available at the upper (cranial) edge of the back plate so that the patient can sit up?

9. Instruct the patient in putting the BORT Generation Osteoporosis Back Brace on.

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BORT GmbH | Postfach 1330 | D-71367 Weinstadt | www.bort.com